


# PHILIPS

## AVENT

### Easy Comfort Single Electric Breast Pump

Includes one 4oz/120ml bottle

 Product Picture

SCF301/05



## Convenience and Comfort for effective expression

Choose between 4 stimulation modes and 4 pumping settings to make every expression session comfortable and personalized. The lightweight, compact design also means expressing on-the-go is a breeze.

### **Lightweight and compact design**

- Easy expressing, at-home or on-the-go
- Battery and electricity operated

### **Multiple hassle-free settings for personalized use**

- Hassle-free, set and forget settings
- 4 stimulation modes x 4 pumping modes for optimal expressing

### **Discreet and hygienic expressing**

- Discreet sound & backlit display for discreet use, day/night

### **Rhythmic Wave Technology**

- Designed to mimic babies' suckling rhythm

### **Designed to gently stimulate milk flow**

- Soft massage cushion designed to help stimulate milk flow

# Highlights

## Personalized pumping

It's important to be comfortable when expressing, and comfort means something different for every mom. With this breast pump you can choose between 4 stimulation modes, and 4 pumping settings. This way, you can make each expressing session efficient and comfortable.

## Hassle-free settings

Moms have enough on their mind without having to recall their ideal setting – so we'll do that for you. Your preferred settings are stored, so at the press of a button you're ready to go. Need a break from pumping? Simply press the pause button, and to resume press play. It's that easy.

## Soft massage cushion

Our breast pumps were designed with your comfort in mind. The soft and velvety texture of the massage cushions feel warm and comfortable against your skin, which can help stimulate milk flow. And the flexible material is designed to mimic your baby's suckling action, which can help milk letdown.

## Lightweight and compact design

The lightweight and compact design means you can easily express at home, or take it with you for expressing on-the-go.

## Dual Power

The product runs on electricity. For extra mobility you can also insert batteries to allow cordless use. You can insert four non-chargeable 1.5V AA batteries or four rechargeable 1.2V AA batteries

## Designed to mimic suckling

When babies breastfeed, they follow a rhythm of suck, swallow, breathe, with a ratio of 1:1:1. This breast pump features special Rhythmic Wave Technology that's designed to mimic a baby's sucking and swallowing pattern, pace and ratio (1:1), to quickly and efficiently generate milk flow. This means your expressing sessions feel natural and comfortable.

## Discreet expressing

Our handy backlit display means you can express discreetly at night, without having to turn on the lights. And the soft sound means less disturbance.

# Specifications

## Design

**Breast pump design:** Compact design  
**Bottle design:** Ergonomic shape, Wide neck

## Country of origin

China

## Material

**Bottle:** BPA free\*, Polypropylene  
**Nipple:** BPA free\*, Silicone  
**Breast pump:** BPA free\* (food contact parts only)

## What is included

**Sealing disc for milk storage:** 1 pcs  
**Breast pump body:** 1 pcs  
**Base unit incl. tubing:** 1 pcs  
**Travel cover:** 1 pcs  
**Natural bottle 4oz:** 1 pcs  
**Soft cushion:** 1 pcs  
**Lanyard:** 1 pcs

## Ease of use

**Breast pump use:** Intuitive assembly, Easy cleaning, Fully compatible range

## Functions

**Rhythmic wave technology:** Mimics baby suckling pattern

## Development stages

**Stage:** 0 – 6 months

