

# PHILIPS

breast pump handle  
including motor unit

ISIS

single electric breast pump

SCF932

 Product Picture

## To adjust the suction rate of the pump

Find matching products on the specifications tab

This control handle adjusts the pumping suction speed and pattern of your Avent breast pump. The pump learns and then repeats your rhythm after only two compressions. The handle is appropriate for single electric breast pumps only.

### **Ease of using Philips spare parts**

- Easy renew your product with original Philips parts

## Highlights

### Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

## Specifications

### Replaceable part

**Fits product types ::** SCF302, SCF312, SCF292, SCF292/01, SCF292/13, SCF312/60, SCF302/11, SCF292/99, SCF302/01, SCF312/13, SCF292/31, SCF302/60, SCF302/13, SCF312/01, SCF292/15, SCF302/19

