

# PHILIPS

## AVENT

Digital baby  
thermometer set

SCH540/01

 Product Picture

## Gentle and accurate

Combining the benefits of a digital thermometer and a pacifier thermometer.

### **Convenient for you and your baby**

- Includes Avent Freeflow Pacifier

### **Fast and reliable digital temperature readings**

- Clear measurement displays and battery status indicators
- Battery low indicator

### **Comfort for your child**

- Soft, flexible tip digital thermometer
- Choice to adapt temperature measurement method to suit baby
- Specially designed pacifier thermometer

# Highlights

## Choice of measurement method

The Philips Avent thermometer accurately and conveniently measures your baby's body temperature orally, axially and rectally.

## Soft, flexible tip digital thermometer

Soft, flexible tip digital thermometer

## Specially designed pacifier thermometer

Enables easy transition from pacifier to temperature measurement.



## Clear measurement displays

Clear measurement displays and battery status indicators

## Battery low indicator

Battery low indicator

## Avent Freeflow Pacifier

Lets your baby get used to the shape.

# Specifications

## Country of origin

Made in China: Thermometers, Yes

## Power

Operating time on battery: >100 hour(s)

Battery type: LR41

Number of batteries: 1 per thermometer

## Accessories

User manual

## Technical specifications

Operating temperature range: 90-109.9 °F

Measurement time: 10secs ±3secs in water (bath), 25secs ±10secs oral & rectal, 60secs ±10secs underarm

Accuracy: 95-107.6 +/-0.2 °F

## Convenience

Automatic shut-off

## Country of origin

Made in England: Soothers, Yes

## Convenience

Flexible tip

Suitable for newborns

Clear measurement display

Battery status indicator

