

PHILIPS

Digital baby
thermometer set

SCH530/86



Easy and reliable temperature measurement

Digital thermometer pacifier for easy and accurate measurement.

Fast and reliable measurement

- Digital thermometer for easy and accurate measurement
- Quick and reliable temperature measurement within seconds

Convenience for you and your baby

- Dimpled shield prevents skin irritation
- Includes a dummy pacifier in the size for a newborn baby
- Designed especially for newborn babies

Baby's safety

- Teat suitable for sterilizing or boiling
- Allows boiling or sterilizing without damaging electronics

Highlights

Digital body thermometer

Fast, reliable and accurate digital thermometer for measuring the baby's temperature.

Measure within seconds

Never too hot, never too cold – the "Just Right" digital thermometer quickly tells you if the milk or food is safe to eat or drink. It takes the guessing out of warming bottles or jars and is great for checking if food is still warm enough after taking a short break in feeding.

Dimpled shield on pacifier

Dimples make sure air flows behind the teat to protect your baby's sensitive skin from irritation.

Silicon teat

The pacifier's specially shaped silicon teat is designed to fit your baby's mouth and will not inhibit natural development. Use it with confidence to pacify your baby.

Dummy pacifier included

You can put the soother regularly in your baby's mouth to stimulate him or her to get used to it. It works like any regular soother.

Detachable electronics

By detaching the sensitive temperature measuring electronics from the pacifier, you can safely sterilize or boil the part your baby comes into contact with.

Sterilize teat

Thanks to the special material, you can safely sterilize the pacifier's teat in a microwave oven or electrical sterilizer or boil it without harming the pacifier.

Specifications

Power

Operating time on battery: 200 hour(s)

Battery type: LR41

Number of batteries: 2

Removable/replaceable

Technical specifications

Operating temperature range: 32–44 °C

Accessories

User manual

Battery

Technical specifications

Accuracy: 36.5–39.0°C:±0.1°C

Measurement time: 10 s

Convenience

Automatic shut-off

Ready indicator

Convenience

Flexible tip

Suitable for newborns

