

PHILIPS

Infrared lamp

HP3614

Relieves muscular pains and common colds

Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth is comfortable and penetrates deep into skin, improving blood circulation thereby giving relief for focused areas of 20x30 cm.

Philips - the light expert

- More than 15 million infrared lamps sold

Highlights

Philips - the light expert

Philips is one of the global leading companies in light, with over 100 years of light experience. Globally Philips already sold more than 15 million of Infrared lamps.

Specifications

Easy to use

Adjustability: 0-40 degrees (backwards)

Technical specifications

Lifetime of lamp: 750 sessions of 10 minutes

Safety

IEC certified: In line with IEC 60335

Easy storage

Cord storage

Technical specifications

Power: 100 W

Cord length: 1.6 m

Voltage: 220/230 V

Frequency: 50 Hz

Insulation: Class II (double isolation)

Type of lamps: R95, 100W

Weight and dimensions

A-box dimensions: 22.4x26.2x58.8
(HxWxD) cm

Product dimensions: 20.0x12.0x16.0
(HxWxD) cm

Product weight: 0.5 kg

F-box weight: 1 kg

No. of F-boxes in A-box: 6 pcs

Qty. on Euro pallet: 288 pcs

A-box weight: 3.1 kg

F-box dimensions: 21.1x12.6x19.2 (HxWxD) cm

Logistic data

Country of origin: Germany

CTV code: 884361401000

