

Philips Sonicare Power  
Flosser 7000  
Oral Irrigator

Quad Stream technology  
Whisper quiet design  
Pulse waves guide as you go



HX3911/41

Removes up to 99.9% plaque in treated areas\*

An effortless, thorough clean

Enjoy an effortless, thorough clean between your teeth and improve your gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth. Easy.

**An effortless clean with advanced technology**

- Easy to fill 600ml reservoir
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Pulse Wave technology guides from tooth to tooth
- Comfort nozzle for sensitive gums
- Whisper quiet for great results without disturbing others
- Quad Stream nozzle for maximum coverage
- Easy magnetic dock with auto shut-off
- Compact nozzle case

**Customize easily to your needs**

- 4 flossing modes, 10 intensities

**PHILIPS**  
**sonicare**

# Highlights

## Covers more area, effortlessly



X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

## Guides from tooth to tooth



Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

## Whisper quiet design




Engineered for a discreet clean, so you get great results without disturbing others.

## 4 modes, 10 intensities




Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean; Deep Clean uses distinct pulses for a more thorough clean; Sensitive mode is extra gentle for sensitive gums; Massage mode soothes and stimulates the gumline. Whenever you need, adjust the intensity for your comfort using the LED light display as your guide.

## F3 Quad Stream nozzle




For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

## F2 Comfort nozzle



Want to take care of sensitive gums? Click on the Comfort nozzle and gently guide the soft rubber tip along your gumline.

## F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

## Magnetic dock w/auto shut-off



When you're finished, simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

## 600ml reservoir



The 600ml reservoir with easy flip lid holds enough water for one cleaning. To refill simply lift out, fill, and click back in.

## Compact nozzle case

Store extra nozzles hygienically in our compact and stylish nozzle case.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

## Performance

- Cleaning: Clean entire mouth in 60-90sec

## Items included

- F1 Standard nozzle: 1
- F3 Quad Stream nozzle: 2

- F2 Comfort nozzle: 1
- Power Flosser: 1
- Nozzle storage case: 1

## Modes

- Clean: For exceptional everyday clean
- Deep Clean+: For an invigorating deep clean
- Intensities: 10
- Massage: Invigorating massage
- Sensitive: Gentle teeth and gum cleaning

## Health benefits

- Gum health: Improves in 2 weeks\*\*
- Plaque removal: Removes up to 99.9% plaque\*, Removes up to 2x more plaque\*\*



Issue date 2023-07-17

Version: 2.2.2

EAN: 87 10103 95331 9

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* Up to 99.9% plaque removal in treated areas in a lab study. Actual in-mouth results may vary.

\*\*With Quad Stream tip on setting 8 with a manual toothbrush.  
Individual results might vary.