




Philips  
Infrared muscular pain  
relief

**200W**  
Focused treatment area



HP3621/02

## Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 200 Watt infrared warmth is very comfortable and penetrates deep into skin, improving blood circulation thereby giving effective relief for focused areas of 20x30 cm.

**PHILIPS**

# Specifications

## Easy to use

- Adjustability: 0-40 degrees (backwards)
- On/off switch
- Soft-touch handgrip: for easy transport and positioning

## Focused treatment

- Treatment area: 20x30 cm
- Infrared halogen lamp: 200 W

## Safety

- IEC certified: In line with IEC 60601 and 60335

## Effective pain relief

- Deep penetrating warmth
- Even warmth distribution: Now, upto 3 x more comfortable warmth

## Technical specifications

- Voltage: 110 or 220 or 230 V
- Cord length: 2 m
- Power: 200 W

- Frequency: 50/60 Hz
- Insulation: Class II (double isolation)
- Lifetime of lamps: 500 hour(s)

## Weight and dimensions

- A-box dimensions: 29x40.4x74.5 (WxHxD) cm
- Product weight: 1.4 kg
- Product dimensions: 20x30.9x17.2 (WxHxD) cm
- F-box weight: 1.8 kg
- No. of F-boxes in A-box: 3 pcs
- A-box weight: 5.4 kg
- Qnt. on Euro pallet: 60 pcs
- F-box dimensions: 28x38.6x24.5 (WxHxD) cm

## Logistic data

- Country of origin: Germany
- CTV code: 884362102000

## Medical appliance

- Proven health benefits: supported by Medical Device Directive certificate



Issue date 2022-09-14

Version: 0.0.0

© 2022 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)