

# PHILIPS

 Product Picture

InfraCare infrared  
lamp

650W

Half body treatment

63.5cm high

HP3641

## Effective pain relief

### Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared warmth is very comfortable and penetrates deep into the skin, stimulates blood circulation and gives effective relief for half-body areas of 60x40 cm.

#### **Easy to use**

- Adjustable angle with rotation

#### **Effective pain relief**

- Deep penetrating Infrared warmth

#### **Half-body treatment**

- 650 Watt infrared halogen lamp

# Highlights

## 650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

### Adjustable angle with rotation

The InfraCare HP 3643 is easy to use for half-body-treatment and for different situations (e.g. lying on a bed or couch) due its maximum flexibility. Next to the height adjustment, the

lamphouse can easily be positioned to optimally treat the right body part, by moving the lamphouse into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

### Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

# Specifications

## Easy to use

**2 handgrips:** for easy transport and positioning

**On/off switch**

**Extendable height:** 45–80 cm

**Easy positionable lamphouse (left/right):** 90

**Easy positionable lamphouse (up/down):** 45

## Half-body treatment

**Infrared halogen lamp:** 650 W

**Treatment area:** 60x40 cm

## Safety

**IEC certified:** Complies to IEC 60601 and 60335

## Effective pain relief

**Deep penetrating warmth**

## Treatment control

**Digital timer:** with auto shut-off

**Dimmer with 3 warmth settings**

**Pause function**

## Technical specifications

**Voltage:** 230 V

**Power:** 650 W

**Cord length:** 300 m

**Frequency:** 50/60 Hz

**Insulation:** Class II (double isolation)

**Lifetime of lamps:** 500 hour(s)

## Weight and dimensions

**Product weight:** 5.5 kg

**Product dimensions:** 29.0x63.5x29.0

(WxHxD) cm

**Qty. on Euro pallet:** 18 pcs

**F-box dimensions (=A-box):** 33.9x80.0x36.5

(WxHxD) cm

**F-box weight (=A-box):** 6.5 kg

## Logistic data

**Country of origin:** Germany

**CTV code:** 884364101000

## Easy to use

**Digital user interface:** for feedback on time and warmth intensity

## Medical appliance

**Medical Device Directive:** 2007/47/EC, MDD 93/42/EEC

