

# PHILIPS

CPAP with Humidifier  
and Heated Tube

DreamStation

Colorful, easy-to-use display

Display and track nightly use

Designed to reduce noise

Bluetooth enabled connectivity



HH1461/00



## Take charge with ease and confidence


With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customize and continue your obstructive sleep apnea (OSA) therapy.

### Performance

- Support for the long term
- Make progress every day
- Small and lightweight
- Performance and appeal


# Highlights

## Bluetooth enabled



Built-in Bluetooth makes connectivity to the DreamMapper self-management tool easy. DreamMapper provides goal-setting tools, helpful videos and ongoing feedback to keep you motivated and engaged.

## Daily Progress Feedback




Each morning you are greeted with Daily Progress Feedback - a summary of how far you've come to motivate you to stick to your therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of the last 30 days of good nights with more than four hours of use.

## Small and lightweight



DreamStation's top-mounted humidifier hose swivel and front-facing display make it easy for you to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

## Performance and appeal



DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

# Specifications

## Humidifier Details

**Tubing compatibility:** Includes heated tubing

## Performance data

**Flex mode:** C-Flex

**Pressure mode:** CPAP

**Ramp time:** 0-45 min (5-min increments)

**Starting ramp pressure:** SmartRamp

## Operating specifications

**Altitude compensation:** Automatic

**Operating pressure:** 4 to 20 cmH2O

## Compatibility

**For use with:** DreamMapper

## Product details

**Bluetooth:** Built-in

