

PHILIPS

Not for sale in the
United States

DreamStation CPAP with Humidifier

Positive airway pressure therapy

Prescription device

Home or hospital use

Not for sale in the United States



HH1460/00



DreamStation CPAP with Humidifier

The Philips Respironics DreamStation systems deliver positive airway pressure therapy for the treatment of Obstructive Sleep Apnea in spontaneously breathing patients weighing over 30 kg (66 lbs). It is for use in the home or hospital/institutional environment.

DreamStation 2

- Patient safety is our top priority

Highlights

Patient safety is our top priority

In the US, Philips Respironics will continue to service sleep and respiratory care devices already with healthcare providers and patients, and supply masks and other accessories, consumables and replacement parts. Until agreed with the FDA, Philips Respironics will

not resume selling new CPAP or BiPAP sleep therapy devices or other respiratory care devices.

Specifications

Performance data

Flex mode: C-Flex

Pressure mode: CPAP

Ramp time: 0-45 min (5-min increments)

Starting ramp pressure: SmartRamp

Operating specifications

Altitude compensation: Automatic

Operating pressure: 4 to 20 cmH2O

Compatibility

For use with: DreamMapper

Product details

Bluetooth: Built-in

