

PHILIPS



Sonic electric
toothbrush

1100 Series

Sonic technology

QuadPacer and SmarTimer

Slim ergonomic design

14-day battery life



HX3641/01

Goodbye manual toothbrush. Hello Sonic technology.

Up to 3x better plaque removal*

Sonic technology combined with our brush action gently removes plaque up to 3x better* than a manual toothbrush.

Designed around you

- 14-day long lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use
- Goodbye adapter. Hello USB.

Optimize your brushing

- Optimize your brushing with SmarTimer and QuadPacer

Easy transition to Sonicare

- Easy-Start to ease the transition

Protects your gums

- Safe and gentle to use

Provides superior clean

- Up to 3x better plaque removal vs a manual toothbrush*
- Our unique technology brings you a powerful yet gentle clean


1100 Series

Sonic electric toothbrush

HX3641/01

Highlights

Advanced Sonicare Technology




Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.** 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

Removes plaque 3x better*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3X better* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

Safe and gentle to use



Sonicare toothbrush provides you with a thorough clean without being harsh on your teeth and gums. Gentle yet powerful sonic vibrations provide an exceptional clean while being gentle on your teeth and gums.

SmartTimer and QuadPacer



2 minute SmartTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

Designed around you



Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Goodbye adapter. Hello USB.

Power adapters may seem small, but added together, their plastics make a huge impact on the environment. By switching to USB charger in our Philips Sonicare 1000-4000 Series alone, we can use up to 189,210 kg less plastic every 4 years. That's powerful, right?

1100 Series

Sonic electric toothbrush

HX3641/01

Specifications

Power

Voltage: DC5V

Technical specifications

Operating time (full to empty): 14 days*****

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Mint green

Service

Warranty: 2-year limited warranty

Ease of use

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Items included

Handle: 1100 Series

Brush heads: 1 C1 SimplyClean standard

Charger: 1

Cleaning performance

Performance: Removes up to 3x more plaque*

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

