

# PHILIPS

Under-the-nose  
mask

Philips Respironics  
Amara View

Full face mask

Under-the-nose cushion

Tube-in-front design

Available in multiple sizes



HH1232/02

## Enjoy the view

Amara View's innovative under-the-nose cushion is designed to prevent red marks, discomfort and irritation on the nose bridge. This tube-in-front, full face mask is also designed to allow users to wear glasses, read, or watch tv at bedtime.

### Wide field of view

- Wide field of view
- Quick and convenient release

### Easy to use

- Easy to use

### Innovative design for comfort

- Innovative design for comfort

## Highlights

### Wide field of view



It's easy to wear glasses, read or watch TV before bedtime without a bulky frame or cushion in front of the face and eyes.

### Innovative design for comfort



Amara View's innovative under-the-nose cushion is designed to prevent red marks, discomfort and irritation on the nose bridge that many patients experience while wearing a traditional full face mask.

### Quick and convenient release



The quick-release tube is designed to aid in easy disassembly as you get up throughout the night.

### Easy to use



Amara View is designed to be easy to use with one modular frame for all cushion sizes.

