

# PHILIPS

BiPAP Pro with Humidifier &  
Heated Tube

DreamStation

Colorful, easy-to-use display

Display and track nightly use

Designed to reduce noise

Bluetooth enabled connectivity

HH1485/02



## Take charge with ease and confidence

With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customize and continue your obstructive sleep apnea (OSA) therapy.

### Performance

- Performance and appeal
- Small and lightweight
- Make progress every day
- Support for the long term

# Highlights

## Performance and appeal

DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

## Small and lightweight


DreamStation's top-mounted hose swivel and front-facing display make it easy for you to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

## Daily Progress Feedback

Each morning you are greeted with Daily Progress Feedback - a summary of how far you've come to motivate you to stick to your therapy. The intuitive, color display shows a

simple trend of their nightly hours of use, followed by a summary of the last 30 days of good nights with more than four hours of use.

## Bluetooth enabled



Built-in Bluetooth makes connectivity to the DreamMapper self-management tool easy. DreamMapper provides goal-setting tools, helpful videos and ongoing feedback to keep you motivated and engaged.

# Specifications

## Compatibility

For use with: DreamMapper

## Product details

Bluetooth: Built-in

## Performance data

Flex mode: C-Flex

Ramp time: 0-45 min (5-min increments)

Starting ramp pressure: SmartRamp

## Operating specifications

Altitude compensation: Automatic

Operating pressure: 4 to 20 cmH<sub>2</sub>O

