

# PHILIPS

Not for sale in the  
United States

DreamStation BiPAP  
autoSV w/ Humidifier

Non-invasive ventilatory support

Prescription device

Home or hospital use

Not for sale in the United States

HH1487/01



## DreamStation BiPAP autoSV with Humidifier

The BiPAP autoSV device is intended to provide non-invasive ventilatory support to treat adult patients (>30 kg/66 lbs) with Obstructive Sleep Apnea and Respiratory Impairment caused by central and/or mixed apneas and periodic breathing. This device may be used in the hospital or home.

### **DreamStation 2**

- Patient safety is our top priority

# Highlights

## Patient safety is our top priority

In the US, Philips Respironics will continue to service sleep and respiratory care devices already with healthcare providers and patients, and supply masks and other accessories, consumables and replacement parts. Until agreed with the FDA, Philips Respironics will

not resume selling new CPAP or BiPAP sleep therapy devices or other respiratory care devices.

# Specifications

## Performance data

**Ramp time:** 0–45 min (5-min increments)

**Starting ramp pressure:** SmartRamp

## Operating specifications

**Altitude compensation:** Automatic

**Operating pressure:** 4 to 20 cmH<sub>2</sub>O

## Compatibility

**For use with:** DreamMapper

## Product details

**Bluetooth:** Built-in

