

# PHILIPS

## AVENT

Double Electric  
breast pump

Breast pumps

Double

Advanced

Corded use

Includes travel bag and pouch

SCF393/62



## Fast and incredibly gentle

Natural motion technology for quick milk flow

Inspired by baby's unique combination of suckling and massage, Natural Motion technology powerfully combines suction and nipple stimulation for quick letdown and milk flow. Finding the right nipple is important. See more information below.

### Personalized and easy to use

- Minimal parts and intuitive setup
- Pause/Play function
- Memory function
- Quiet motor for a discreet experience anywhere
- Breast pump bag and pouch

### Gentle and comfortable

- Soft and adaptive silicone cushion
- Express without leaning forward
- Personalized experience, 8 + 16 setting levels

### Quick expressing, more milk in less time\*

- Natural motion technology for a quick milk flow\*

### Finding the right nipple is important

- Finding the right nipple is important

## Specifications

### Material

**Bottle:** BPA free\*, Polypropylene

**Nipple:** BPA free\*, Silicone

**Breast pump:** BPA-free\* (food contact parts only)

### Ease of use

**Breast pump use:** Easy cleaning, Intuitive assembly, memory function

### Functions

**No leaning forward:** Sit in a comfortable position

**Settings:** 16 expression levels, 8 stimulation levels

**Soft massage flange:** Gentle stimulation

### Development stages

**Stage:** 0 - 6 months



\* Based on milk flow initiation time (time to Milk Ejection Reflex - MER) results from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018).

\* 1) Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 2019, (109 participants, Israel); (2) Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

\* Journal of Obstetric, Gynecologic & Neonatal Nursing, 1993, (20 Caucasian participants, USA); (3) Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

\* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011.