

Philips Philips Avent
Silicone tube

Breast pump

CP2061/01



Connects parts of your electric breast pump

Check for compatibility below

The silicone tube runs between the breast pump and motor unit. One tube included per pack. For use with electric breast pumps.

PHILIPS

Silicone tube
Breast pump

CP2061/01

Specifications

Replaceable parts

• Fit product types: SCF332, SCF334, SCF391,

SCF393, SCF394, SCF395, SCF396, SCF397,
SCF398



Issue date 2024-04-13

Version: 4.4.1

© 2024 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com