

# PHILIPS

## AVENT

### Essential

Single Electric  
Breast Pump

SCF323/11



## Express effortlessly

The Philips Avent electric breast pump keeps milk flow at an optimum and gently adapts to your size and shape. Quick, comfortable, and convenient for each new mother.

### Personalized


- Personalised experience, 4 + 4 setting levels
- Minimal parts and intuitive setup
- Memory function
- Small and Compact

### Gentle and comfortable

- Soft and adaptive silicone cushion

# Highlights

## Personalised experience



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 4 stimulation and 4 expression levels for a personalised experience.

## Soft adaptive silicone cushion



One size fits all. The silicone cushion gently flexes and adapts to fit your nipple.

## Easy to clean and set up



Our closed expression system means milk stays out of the tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

## Memory setting



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

## Small and Compact



Small and Compact motor unit

# Specifications

## Power

USB-C connection

Main voltage: 5V and  $\geq 1$  A

Supply unit is not included. Not chargeable, corded use only.

## Material

0% BPA\*

## What is included

Motor unit

Silicone tube

USB-C cable

Natural bottle (125 ml)

Expression kit

Screw-on lid

## Ease of use

**Breast pump use:** Easy cleaning & assembly, memory function, portable design

## Functions

**Settings:** 4 expression levels, 4 massage levels

