



Philips  
Knee massager

**Ultra-thin body**  
3 levels of vibration  
Detachable controller

PPM5301BU

## Enjoy life with healthy joints

Protect knees and care for your joints

The knee massager has hot compress and vibration massage functions, which can effectively promote blood circulation, relieve soreness and discomfort, and care for the knee joints. Detachable controller, light and portable to carry.

### **Combined massage to relieve discomfort**

- Vibration mode, professional relaxation massage
- Hot compress of 40-60°C, 5 adjustable levels.

### **Ergonomic design, fit and comfortable**

- 3D surround ergonomic design for better protect knees

### **Detachable controller, light and portable**

- The controller is detachable, no burden to carry around

# PHILIPS

Knee massager  
Ultra-thin body 3 levels of vibration, Detachable controller

PPM5301BU/90

# Specifications

## General specifications

- Net weight: 280g
- Color: Dark blue

## Main parameters

- Temperature of heating: 40-60 °C
- Rated voltage: 5V

- Battery: 2000mAh
- Charging port: Type-C

## Country of origin

- Made in: China

## Service

- Warranty period: 2 years

# Highlights

## Vibration mode

3 levels of adjustable vibration: level 1 to wake up muscles, level 2 to relax muscles, and level 3 to deep massage. You can choose suitable gears by yourself to avoid damage and effectively remove fatigue.

## Hot compress function

The temperature of 5-level hot compress can be adjusted to promote vasodilation, strengthen local blood circulation, effectively relieve muscle tension, reduce local pain. You will have a full relaxation and fatigue relief of the day.

## 3D wraparound design

Accurately massage the important acupoints of the knee, and cover more massage area. Using skin-friendly fabrics to enhance the experience, promote blood circulation, and relieve discomfort.

## Detachable controller

The main unit is detachable, which is more convenient to carry and clean. The host has a built-in battery with a capacity of 1800mAh, which can be used for 90 minutes when fully charged.



Issue date 2023-04-17

Version: 2.2.1

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)