

Philips
Sports massage gun

Professional grade materials
15kg strong power output
Real-time pressure feedback

PPM7521BK

Strong from inside out

15kg professional thrust at very low noise

The ultra -silent and all-metal power system provides up to 15kg of strong thrust. With 12 high amplitude and 8 professional grade massage heads, including a skin friendly titanium head, this massage gun is idea for deep tissue treatment.

Powerful massage performance

- Durable aluminum body
- Strong thrust output up to 15kg
- Human-computer interaction make massage easier

Easy to use

- Non-slip design handle, make it easy to grasp
- High intensity plastic hardshell case for full protection.

Professional massage components

- 8 different massage heads meet all your needs
- Enhance the massage amplitude and comfort level
- Deliver strong penetrating waves while avoid skin allergies

PHILIPS

Sports massage gun
Professional grade materials 15kg strong power output, Real-time pressure feedback

Specifications

General specifications

- Net weight: 980g
- Product Dimension (LxWxH): 78x142x48 mm
- Main materials: Aluminum alloy, Silicone
- Color: Titanium gray

Main parameters

- Working voltage: 22.2V
- Massage amplitude: 12mm
- Massage frequency: 1600Hz-3200Hz

- Maximum thrust: 15kg
- Battery: 2600mAhx6
- Charging port: DC 26V 1A

Country of origin

- Made in: China

Service

- Warranty period: 2 years

PPM7521BK/79

Highlights

Strong output

SuperDynamic all-metal power system, providing 15kg professional thrust. The massage waveform is more stable and effective.

Aluminum body

The massage gun is made of aviation aluminum alloy, which helps to reduce the weight while ensuring stable structural strength. Robust and durable.

Human-computer interaction

Intelligent strength feedback display, helps to control massage pressure and intensity. For an optimal massage, while avoiding excessive pressure.

Professional accessories

Each massage head is carefully designed taking into account the needs of the various body parts, based on the principles of physiotherapy.

Professional titanium head

The professional titanium head is popular with physiotherapists. The metal massage head with higher hardness enables the massage gun to reach deep into the muscles, target the deep to tissue layers that ordinary massage heads cannot reach and completely soothe muscle groups. Titanium is the most friendly metal to the body and does not cause discomfort to people with skin allergies.

High elastic rubber ball head

The high elasticity of the rubber complements the vibration frequency of the massage gun, optimizing muscle tissue penetration for the best possible massage effect. The soft rubber surface is easy to clean supporting bacteria and harmful substances removal after long-term use. Protecting the skin's health.

Non-slip design, easy grasp

The 45mm slim handle covered with non-slip silicone, supports a convenient everyday massage.

Anti-shock portable case

The small and portable case features a high-strength plastic shell with anti-drop and anti-crush properties for providing full protection for the machine.



Issue date 2022-12-13

Version: 1.0.1

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com