

# PHILIPS

## sonicare

Sonic electric  
toothbrush

Essence

1 mode

1 brush head

HX5361/29



## Better plaque removal

### The ultimate Sonicare toothbrush

The unique dynamic cleaning action of Sonicare's electric toothbrush, gently and effectively reaches deep between teeth and along the gumline.

#### **Proven to improve oral health**

- Clinically proven safe and gentle

#### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time

#### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

#### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle




Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## SmartTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

# Specifications

## Design and finishing

Color: White

## Power

Voltage: 110-220 V

## Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: NiMH

## Service

Warranty: 2-year limited warranty

## Ease of use

Brush head system: Easy screw-on brush heads

Battery indicator: Light shows battery status

Handle: Ergonomic design, Rubber grip for easy handling

Brushing time: Up to 2 weeks

## Items included

Charger: 1

Brush heads: 1 e-Series standard

Handles: 1 Essence

Travel case: 1

## Cleaning performance

Performance: Up to 2X more plaque removal\*

Health benefits: Helps improve gum health

Speed: Up to 31000 brush movement/min

Whitening benefits: Helps naturally remove stains

Timer: SmartTimer

## Modes

Clean: For exceptional everyday clean

