

PHILIPS

Power Toothbrush
Special Edition

9000 Series

Connected app tracks progress

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities



HX9911/90



An effortlessly white smile

Whiter teeth in just 3 days*

Start every day in style with advanced Sonicare technology for whiter teeth and oral health. Gentle on gums with a connected app for coaching and tracking. Choose from 4 modes and 3 intensities for a personalized brushing experience.

Whiter teeth & improved oral health

- Hard on plaque, gentle on your gums

Adapts to your needs

- 4 modes & 3 intensities for personalized brushing experience

Adapts to your needs

- Pressure sensor lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you

Whiter teeth & improved oral health

- Remove surface stains for whiter teeth in just 3 days*
- Advanced Sonicare technology for a powerful yet gentle clean

Highlights

Whiter teeth in just 3 days*



Boost your smile's brightness with advanced Sonicare technology. Polishing bristles remove up to 100% more stain than a manual toothbrush for visibly whiter teeth in just 3 days*.

Gentle on your gums




The C3 Premium Plaque Control brush head is designed to give you your deepest clean yet. Soft, flexible sides and bristles perfectly curve around the contours of your teeth, giving you 4x more surface contact and helping target hard-to-reach spots***

Advanced Sonicare technology



Powerful sonic vibrations move 62,000 times every minute to drive plaque-fighting bubbles between your teeth and along your gum line. Removes up to 10x more plaque than a manual toothbrush.

4 modes / 3 intensities



Choose from 4 brushing modes: Clean for exceptional everyday cleaning, White+ to brighten your smile, Gum Health for a gentle yet effective clean, and Deep Clean+ for an invigorating deep clean. Three intensity settings from Low to High let you boost performance or care for sensitive gums.

Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work.

BrushSync automatic pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. When you are using the C3 Premium Plaque Control brush head, your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the Clean mode to help whiten your teeth.

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 14 days***

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Silk Pink to White Gradient

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Brush heads: 1 C3 Premium Plaque Control

Travel case: 1

Handle: 1 DiamondClean 9000

Charger base and puck: 1

Cleaning performance

Performance: Removes up to 10x more plaque*

Timer: BrushPacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

3 intensities: High, Medium, Low

Clean: For exceptional everyday clean

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

White+: To remove surface stains

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads



* in White+ Mode using a leading whitening toothpaste

** Vs a manual toothbrush, actual results may vary

*** as compared to DiamondClean brush head