

**PHILIPS**

**AVENT**



Large massage  
cushion

Fits all comfort breast pump



SCF167/01

## Right size, maximum comfort

### Breast pump with wider funnel

This soft massage cushion has a wider funnel, specially designed for mums with larger nipples. It can be used with Philips Avent Comfort breast pump range. See below images to help you decide if you need the larger breast pump cushion.

#### **Comfortable pumping**

- Wider funnel for comfortable fit
- Fits all Natural range breast pumps

#### **Check your fit for optimal comfort**

- Right fit: The nipple sits comfortably in the cushion funnel
- Choose large size: Nipple touches the walls of the funnel

# Highlights

**Large breast pump cushion**

You can find out whether you need a large breast pump cushion by checking whether your nipple fits comfortably in the cushion funnel. If your nipple touches or rubs against the walls of the cushion funnel, then you should choose the large breast pump cushion to ensure a more comfortable expressing experience.

**Correct fit**

If the nipple sits comfortably in the cushion funnel, you are using the correct size massage cushion.

**Too small fit**

If your nipple touches or rubs against the wall of the cushion funnel, then you should choose the large breast pump cushion with wider funnel to ensure a more comfortable pumping experience.

**Fits all Natural breast pumps**

The large size soft massage cushion is compatible with all Philips Avent Comfort breast pumps in the Natural range.

# Specifications

**Cushion**

**Material:** Silicone, BPA free\*

**Country of origin**

England

**Material**

BPA free\*

**What is included**

Large size cushion: 1 pcs

**Compatibility**

Compatible with Natural range

**Ease of use**

**Cushion:** Easy to clean, Dishwasher proof

**Functions**

**Soft massage cushion:** Gentle stimulation

**Development stages**

**Stage:** Pregnancy, 0 - 6 months

**Replaceable part**

**Fits product types ::** SCF330, SCF332, SCF334, SCF330, SCF332

