

# PHILIPS

## sonicare



Sonic electric  
toothbrush

Optimal Clean

Built-in pressure sensor

2 cleaning modes

1 x BrushSync feature  
with sanitizer



HX6829/76

## Whiter teeth. Gently does it.

Whitens teeth in just 1 week.

Feel the difference of a gentle clean with our Pressure Sensor while whitening your teeth in 1 week

### Proven to improve oral health

- Up to 7x more plaque removal than a manual toothbrush
- Safe & gentle on sensitive areas, orthodontics & dental work
- Whitens teeth in just one week
- Clean and white mode

### Designed around you

- Makes traveling easier for you

### Innovative technology

- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads

## Specifications

### Power

Voltage: 110-220 V

### Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: Lithium ION

### Design and finishing

Color: White and Light Blue

### Service

Warranty: 2-year limited warranty

### Ease of use

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Brushing time: Up to 2 weeks

Brush head recognition: Syncs with optimal mode

Replacement reminder: To always ensure best results, reminder icon lights up

### Items included

Handle: 2 ProtectiveClean

Brush heads: 2 W2 Optimal White standard, 2 C2 Optimal Plaque Control

Travel case: 1 hard travel case

Charger: 2

UV sanitizer: 1

### Cleaning performance

Health benefits: Helps reduce cavities

Performance: Removes up to 7x more plaque\*

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Whitening benefits: Whitens teeth in 1 week

Pressure feedback: Vibrates handle to alert user

### Modes

Clean: For exceptional everyday clean

White: Removes surface stains

### Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

