

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

EasyClean

1 mode

2 brush heads

HX6582/51



Healthier gums in 2 weeks

The switch from a manual toothbrush made easy

Philips Sonicare EasyClean combines our patented Sonicare toothbrush technology and the ProResults gum health brush head to remove up to 6x more plaque*. It offers an effective solution to prevent and reduce the early signs of gum disease.

Designed around you

- EasyStart
- Recharge indicator

Guides you to follow dental recommendations

- Quadpacer
- Smartimer

Provides a superior clean

- Designed for optimal gum health improvement
- Patented Sonicare toothbrush technology
- Exclusive bonus items only available to dental professionals

Highlights

Bonus items not sold at retail



Now you can offer an exclusive Philips Sonicare EasyClean bonus item your patients won't find anywhere else, for improved gum care. Designed around patient needs for dental offices like yours, the package includes: 1 extra ProResults gum health brush head.

A gentler experience




This Philips Sonicare toothbrush head features densely-packed, high-quality tufts to effectively remove plaque and improve gum health. Specially curved power tufts to easily reach teeth at the back of the mouth. Also features gentle trim profile to cushion teeth for an extra gentle brushing experience.

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

EasyStart



Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.

Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

Smartimer



Two-minute timer helps ensure dental professional recommended brushing time.

Recharge indicator

Light blinks to let user know the battery needs recharging.

Specifications

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Glacier Green

Service

Warranty: 2-year limited warranty

Ease of use

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Items included

Charger: 1

Brush heads: 2 ProResults gum health

Handles: 1 EasyClean

Cleaning performance

Performance: Removes up to 6X more plaque*

Health benefits: Improve gum health in 2 weeks*

Speed: Up to 62000 brush movement/min

Whitening benefits: Helps maintain white teeth

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

