




Philips Sonicare  
AirFloss Pro/Ultra -  
Interdental cleaner

w/ 1 nozzle & 1 brush head  
w/ Gum health toothbrush



HX8392/45

## Healthier gums in 2 weeks, guaranteed\*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra\*\*\*\* is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.\*\*

### Improves oral health

- Clinically proven as effective as floss for gum health\*\*
- Up to 99.9% plaque removal\*\*\*
- Helps prevent cavities between teeth

### Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings


### The easiest way to effectively clean between teeth

- An easy way to start a healthy habit
- Total daily confidence of a fresh and clean mouth
- Takes just 60 seconds: Point, Press, Clean!

**PHILIPS**  
**sonicare**


# Highlights

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*

## Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

## Fresh & clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

## Helps prevent cavities




By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

## High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

# Specifications

## Design and finishing

- Color: Gray

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Consumer Trade Item

- Height: 10.3
- Width: 11.2
- Length: 22.2
- Net Weight: .52
- Gross Weight: .669
- GTIN: 08710103844914
- Country of origin: MX
- Harmonized system code: 850980

## Outer Carton

- Height: 23.8
- Width: 31.7
- Length: 23.2
- Gross Weight: 4.334
- GTIN: 18710103844911

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Battery Life: 1-3 weeks\*\*\*\*
- Nozzle attachment: Easily snaps on and off

## Items included

- Brush heads: 1 gum health brush head
- Charger: 1
- Rechargeable toothbrush handle: 3 Series gum health
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 1

## Cleaning performance

- Performance: Removes up to 99.9% plaque\*\*
- Health benefits: Cleans entire mouth in 60 sec, Improve gum health in 2 weeks\*
- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
- For best results: Change nozzle every 6 months



Issue date 2023-01-14

Version: 2.2.1

EAN: 08 71010 38449 14

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.  
\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary  
\* or your money back  
\*\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.  
\*\*\*\*\* depending on burst setting used