

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

For Kids

2 modes

2 brush heads

10 stickers



HX6382/07



Get the most from each brushing

A healthy smile for life

Brushing techniques take time to develop. That's why most traditional toothbrushes don't make the cut. But Philips Sonicare For Kids removes up to 75% more plaque in hard-to-reach places than a manual brush.*

Designed around you

- Two brush head sizes
- Ergonomic, easy-grip handle
- Eight faceplate stickers
- Two kid-friendly power modes
- Rubber overmold
- Two brush head sizes

Guides you to follow dental recommendations

- KidTimer

Provides a superior clean

- Patented Sonicare toothbrush technology

Highlights

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

KidTimer



Gradually increases brushing time over 90 days until reaching the dentist-recommended two minutes to train kids to brush longer.

Two brush head sizes



Provide the perfect option for every mouth. Standard and compact sizes are comfortable fit for kids ages 4 and older.

Brush heads



Provide the perfect option for every mouth. Standard and Compact sizes are comfortable fit for kids ages 3 and older.

Power modes



Give them just the right amount of power — low mode for younger children and high mode for older kids.

Faceplate stickers

Make it fun and customizable. Kids can make each brush uniquely their own.

Rubber overmold



Provides added protection for growing teeth and a safer, gentler brushing experience.

Ergonomic handle



For small hands to maneuver more independently and comfortably.

Specifications

Design and finishing

Color: Aqua

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 3 week battery life

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Display: Illuminated display

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design, Rubber grip for easy handling

Brushing time: Up to 3 weeks**

Items included

Handles: 1 Sonicare for Kids

Charger: 1

Brush heads: 1 Sonicare for Kids standard, 1 Sonicare for Kids compact

Stickers: 10 customization stickers

Cleaning performance

Performance: 75% more effective**

Health benefits: For healthy oral care habits

Speed: Up to 62000 brush movement/min

Timer: KidTimer and Quadpacer

Modes

Power modes: 2

