

Philips Sonicare AirFloss
Sonicare AirFloss

Rechargeable
w/ 1 nozzle



HX8111/12

An easier way to floss

Gently and effectively improves oral health

If flossing isn't a regular part of your oral care routine, you're not reaching the bacteria that can build up between teeth causing infection. Sonicare AirFloss helps you to get a deep clean everyday by taking the hassle out of flossing.

Achieve a 30-second clean

- Simple one-button use; as easy as a mouse-click

Proven improvement of oral health

- Up to 99.9% plaque removal

Simple: point and press

- Easy-to-fill reservoir
- Guidance tip for easy placement
- Microburst technology cleans at the touch of a button
- Slim, angled nozzle makes it easy to clean even the hard-to-reach areas
- Slim, easy to hold handle
- Fits standard Sonicare chargers
- Charging after two weeks

PHILIPS
sonicare


Highlights

Microburst technology



Microburst technology applies a quick burst of pressurised air and micro-water droplets to clean deep between teeth where a toothbrush can't reach

Slim, angled nozzle



Slim, angled nozzle makes it easy to clean even the hard-to-reach areas

Slim, easy to hold handle



Easy to manoeuvre around the mouth.

Charging after two weeks



The AirFloss lasts two weeks between charging.

Guidance tip



Simply glide the guidance tip along the gum-line until you feel it settle between the teeth. Provides easy positioning even on those hard-to-reach back teeth.

Fits Sonicare chargers



You only need one charger to charge your Sonicare products.

Simple one-button use




Simple one-button use for easy cleaning in 60 seconds.

Easy-to-fill reservoir



Easy-to-fill reservoir uses only a tiny amount of water for less mess. Fill with mouthwash for enhanced oral health benefits.

Up to 99.9% Plaque Removal



AirFloss uses microdroplets of air and water or mouthwash to remove up to 99% more plaque between teeth*

Specifications

Power

- Voltage: Multi-voltage charger

Technical specifications

- Operating time (full to empty): 2 weeks or 14 full-mouth uses
- Battery: Rechargeable
- Operating time: Replace nozzle every 6 months
- Charging time: 24 hour(s)

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Light in power button blinks to indicate when to recharge



Issue date 2023-01-14

Version: 2.0.1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Compared with manual tooth brushing alone