

Philips Zoom Ultimate
Protocol
In-office whitening
treatment



DIS998/11

A proven whitening solution for deep stains

Your dental professional now has a complete whitening strategy to effectively treat deep staining. The new Philips Zoom Ultimate Protocol combines Philips Zoom WhiteSpeed chairside whitening with Philips Zoom NiteWhite take-home treatment.

Professional-quality results

- Case study findingsCase study
- Superior results for the toughest stains
- What to expect

PHILIPS ZOOM!

Specifications

Performance

- Customizable
- Dual-barrel syringe: No refrigeration required
- Sensitivity: Advanced sensitivity relief
- Maximum whitening: Improved luster

Items included

- Dual-barrel syringes: 2 NiteWhite
- Tray case: Holds custom whitening tray
- Chairside kits: 2
- Light guide: 1
- Relief ACP[®] Oral Care Gel: 1
- Travel bag: 1
- Safety goggles: 1
- Whitening LED accelerator: 1

Proven safe

- Sensitivity management: Adjustable light intensity
- Safe: Little to no sensitivity

Settings

- Intensity: Low, medium, high

Formula components

- ACP: Amorphous calcium phosphate
- Potassium nitrate
- Fluoride
- Carbamide peroxide: 16%

Highlights

For the toughest stains

Philips Zoom Ultimate Protocol consists of four 15-minute Philips Zoom WhiteSpeed sessions, followed by two to four take-home treatments of Philips Zoom NiteWhite Maximum white to achieve your desired whitening results.

Case study findings

In a Philips Zoom Ultimate Protocol case study, patients saw an average improvement of 11 shades and several patients saw results of 14 shades whiter.*

What to expect

At your initial consultation, your dental professional will verify that you're a candidate for whitening, assess your shade and take impressions for your custom fabricated trays. Your in-office procedure will begin with prep and the first of four Philips Zoom WhiteSpeed in-office whitening gel applications. The WhiteSpeed lamp accelerates whitening so that each of the four sessions last only 15 minutes. Post treatment your dental professional will provide custom trays along with the Zoom DayWhite or NiteWhite best for you so you can achieve a bright, white smile.



Issue date 2023-01-14

Version: 4.4.3

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Case study results may not be representative of actual results