




Philips Sonicare  
AirFloss Pro/Ultra -  
Interdental cleaner

w/ ProtectiveClean,  
1 nozzle  
1 brush head



HX8424/17

## Healthier gums in 2 weeks, guaranteed\*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra\*\*\*\* is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is as effective as floss for gum health.\*\*

### Improves oral health

- Up to 99.9% plaque removal\*\*\*
- Helps prevent cavities between teeth

### Optimize your brushing

- Up to 7x more plaque removal than a manual toothbrush
- As effective as floss for gum health\*\*

### Superior technology

- High performance nozzle
- Triple burst customizable settings
- Connects smart brush handle and smart brush heads
- Air and micro-droplet technology

### The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- Total daily confidence of a fresh and clean mouth
- An easy way to start a healthy habit

**PHILIPS**  
**sonicare**

# Highlights

## Say goodbye to plaque



Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush. And the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

## BrushSync technology




A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*

## Micro-droplet technology




The results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Better gum health



Philips Sonicare AirFloss Pro/Ultra can help improve gum health as much as floss\*\* in as little as two weeks.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

## Fresh & clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

# Specifications

## Design and finishing

- Color: Pink, White with grey accents

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Battery Life: 1-3 weeks\*\*\*\*
- Nozzle attachment: Easily snaps on and off

## Items included

- Handle: 1 ProtectiveClean
- Charger: 1
- Brush heads: 1 C2 Optimal Plaque Control
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 1

## Cleaning performance

- Health benefits: Improve gum health in 2 weeks\*, Cleans entire mouth in 60 sec
- Performance: Removes up to 99.9% plaque\*\*
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
- For best results: Change nozzle every 6 months
- Timer: Quadpacer and SmarTimer

## Modes

- Clean: For exceptional everyday clean
- 2 intensities: High, Low



Issue date 2023-01-14

Version: 2.2.1

EAN: 08 71010 39073 81

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* or your money back

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary

\*\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

\*\*\*\*\* depending on burst setting used

\* Removes up to 7x more plaque vs. a manual toothbrush