

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

2 modes

1 brush head

HX6711/02



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush HX6711/02 is proven to remove everyday stains in just two weeks, with the use of the Clean & White Mode.

Proven to improve oral health

- Philips Sonicare toothbrush helps whiten teeth

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic sonic action drives fluid between teeth
- Removes up to 2x more plaque than a manual toothbrush

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains


Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Helps whiten teeth




This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth



Specifications

Design and finishing

Colour: Glacier Green

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 3 weeks**

Battery: Rechargeable

Battery type: Lithium ION

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design, Rubber grip for easy handling

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks**

Items included

Brush heads: 1 ProResults standard

Charger: 1

Handles: 1 HealthyWhite

Cleaning performance

Performance: Removes up to 90% of stains*

Speed: Up to 62,000 brush movements/min

Health benefits: Helps improve gum health

Whitening benefits: Whitens teeth by up to 2 shades

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

Clean and White: 2 min Clean + 30 sec White

