

# PHILIPS

## sonicare

Sonic electric  
toothbrush

Essence

1 mode

2 brush heads

HX5610/04



## Removes up to 2X more plaque\*

\* vs a manual toothbrush

The unique dynamic cleaning action of Sonicare's electric toothbrush, gently and effectively reaches deep between teeth and along the gumline.

### Proven to improve oral health

- Clinically proven safe and gentle
- Philips Sonicare toothbrush helps whiten teeth
- Removes up to 2x more plaque than a manual toothbrush

### Designed around you

- Soft grip helps provide better hand control & added comfort.
- Easy-start to acclimate to the Sonicare experience

### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time

### Provides a superior clean

- Angled brush head neck for better reach to back teeth
- Dynamic cleaning action drives fluid between teeth

# Highlights

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## SmarTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Soft Grip




The soft grip helps you better control your Sonicare toothbrush and adds comfort.

## Patented Sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Easy-start



Gently increases power over first 12 uses, to ease into the Sonicare experience

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: NiMH

## Design and finishing

Color: White with blue soft grip

## Service

Warranty: 2-year limited warranty

## Ease of use

Handle: Ergonomic design, Rubber grip for easy handling

Brush head system: Easy screw-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

## Items included

Travel case: 1

Brush heads: 2 e-Series standard

Handles: 2 Essence

Charger: 1

## Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 31000 brush movement/min

Performance: Up to 2X more plaque removal\*

Whitening benefits: Helps naturally remove stains

Timer: SmarTimer

## Modes

Clean: For exceptional everyday clean

