

# PHILIPS

## sonicare

Sonic electric  
toothbrush

Essence

1 mode

1 brush head

HX5630/50



## Removes up to 2X more plaque\*

\* vs a manual toothbrush

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

### Proven to improve oral health

- Clinically proven safe and gentle
- Philips Sonicare toothbrush helps whiten teeth
- Removes up to 2x more plaque than a manual toothbrush

### Designed around you

- Easy-start program to acclimate to the Sonicare experience
- Soft grip helps provide better hand control & added comfort.

### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time

### Provides a superior clean

- Angled brush head neck for better reach to back teeth
- Dynamic cleaning action drives fluid between teeth

# Highlights

## Patented Sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Smarter




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Removes up to 2x more plaque




Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

## Easy-start program




Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Soft Grip



The soft grip helps you better control your Sonicare toothbrush and adds comfort.

## Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.



# Specifications

## Design and finishing

**Color:** White with red soft grip

## Power

**Voltage:** 110-220 V

## Technical specifications

**Operating time (full to empty):** Up to 2 weeks

**Battery:** Rechargeable

**Battery type:** NiMH

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy screw-on brush heads

**Handle:** Ergonomic design, Rubber grip for easy handling

**Brushing time:** Up to 2 weeks

## Items included

**Handles:** 1 Essence

**Charger:** 1

**Travel case:** 1

**Brush heads:** 1 e-Series standard

## Cleaning performance

**Performance:** Up to 2X more plaque removal\*

**Health benefits:** Helps improve gum health

**Speed:** Up to 31000 brush movement/min

**Timer:** SmarTimer

**Whitening benefits:** Helps naturally remove stains

## Modes

**Clean:** For exceptional everyday clean

