

PHILIPS

sonicare

Sonic electric
toothbrush

CleanCare

1 mode

1 brush head

HX5350/02



Electric toothbrush for better plaque removal

The ultimate toothbrush

The unique dynamic cleaning action of Philips Sonicare electric toothbrush HX5350/02 gently and effectively reaches deep between teeth and along the gum line.

Proven to improve oral health

- Clinically proven safe and gentle
- Philips Sonicare toothbrush helps whiten teeth

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- Sonicare dynamic sonic action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic sonic action drives fluid between teeth

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience

Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Angled brush head neck



The uniquely angled neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle




Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Helps whiten teeth



This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Patented Sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.



Specifications

Design and finishing

Colour: White

Power

Voltage: 110-120 V, 50-60 Hz, 220-240 V, 50-60 Hz

Technical specifications

Battery: Rechargeable

Charging time: 24 hour(s)

Power consumption: 0.9/1.1 W

Service

Warranty: 2-year limited warranty

Ease of use

Brush head system: Easy screw-on brush heads

Battery indicator: Light in power button blinks to indicate when to recharge

Brushing time: 28 2-minute brushings OR 2 weeks

Items included

Charger: 1 charger

Charger base: 1

Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 31,000 brush movements/min

Brush heads: 1 e Series standard brush head

Performance: Better plaque removal than a manual toothbrush

Whitening benefits: Helps remove tea, coffee and other stains

Timer: Smartimer

