

# PHILIPS

## AVENT

Single electronic  
breast pump

Includes 4oz bottle

SCF312/01



## Designed for Comfort

Optimum comfort breast pump with massage cushions

Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. Our Philips Avent breast pump SCF312/01 has been designed to increase your comfort when you are pumping.

### Designed for comfort

- Remembers your personal rhythm
- Patented soft 5-petal massage cushion

### Clinically proven effective

- Gentle yet effective results

### Other great features

- Unique system for easy milk storage
- Also use as a manual breast pump

# Highlights

## Remembers your personal rhythm



Only Philips Avent's electric breast pumps have a unique memory feature that learns and continues your personal pumping rhythm for maximum comfort

## Gentle yet effective results



The gentle-draw vacuum of the Philips Avent breast pump mimics a baby's natural suckling action for comfortable efficacy

## Patented soft massage cushion



The unique active massage cushion is intended to help stimulate natural let-down

## Unique milk storage system



Make your life easier by expressing directly into any of our wide range of Philips Avent baby bottles and breast milk storage containers for the fridge or freezer

## Also u as a manual breast pump



Can also be used as a manual breast pump when that's more convenient

# Specifications

## Power

Voltage: 100 - 240 V

## Country of origin

England

## Material

BPA free\*

## What is included

Breast milk container (125 ml/ 4 oz): 1 pcs

Battery pack: 1 pcs

Bottle stand/funnel cover: 1 pcs

Electronic Breast Pump: 1

Extra soft newborn flow nipple: 1 pcs

Disposable breast pad: 1 pcs

Newborn nipple travel pack: 1 pcs

Spare parts: 2 pcs

Manual pump parts for when power is not available: 1

## Development stages

Stage: Pregnancy, 0 - 6 months

