

PHILIPS

Electric breast pump



SCF396



Inspired by baby, made for mums everywhere

Enter a new era of expressing with a perfect balance of suction and nipple stimulation inspired by the natural way babies drink. The Philips Avent electric breast pump keeps milk flow at an optimum and gently adapts to your nipple size and shape

Personalised and easy to use

- Minimal parts and intuitive setup
- Insulation Pouch included
- Start/Pause function
- Memory function
- Quiet motor for a discreet experience anywhere
- Personalised experience, 8 + 16 setting levels

Gentle and comfortable

- Soft and adaptive silicone cushion
- Express without leaning forward

Quick expressing, more milk in less time*


- Natural motion technology for a quick milk flow*

Rechargeable battery and timer display

- Rechargeable battery for pumping on the go
- Display to keep track of time


Highlights

Quick milk flow*



Express more milk in less time* with a cushion that stimulates the breast to express milk just like baby. It seamlessly adjusts from stimulation mode to expression mode and applies just the right amount of teat stimulation and suction for maximum milk flow. Based on milk flow initiation time (time to Milk Ejection Reflex — MER) results.*

Adapts to your nipple



Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your nipple. It fits 99.98% of nipple sizes* (up to 30 mm).

Sit comfortably



Feel relaxed as you express thanks to a design that lets you sit upright rather than having to lean forward. No need to lean forward, a clinically proven comfortable expressing position*. Based on product clinical trial results with 20 participants (2019); 90% of participants find the expressing position

comfortable (single electric); 95% of participants find the expressing position comfortable (double electric).

Personalise to your needs



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 8 stimulation and 16 expression levels for a personalised experience.

Express anywhere discreetly



Express whenever and wherever you like thanks to our quiet motor with a slimline design.

Use easily on the go



Express on the go without having to worry about a power supply. Our rechargeable battery comes with a micro-USB adapter and lasts up to 3 sessions on a single full charge.

Easy to clean and set up



Our closed expression system means milk stays out of the tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Save your favourite settings



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Track your pumping time

Track how long you express on each breast and your total pumping time with the digital timer display.

Pause whenever you like



If you need to adjust or take a break, the start/pause button is always at your fingertips.

Specifications

Power

Mains voltage: 100V - 240 V

What is included

Expression kit: 1 pcs

Micro-USB adapter: 1 pcs

Motor unit (Rechargeable): 1 pcs

Sealing disc: 1 pcs

Milk bags: 5 pcs

Travel Pouch: 1 pcs

