



Philips Avent
Cap

Silicone Diaphragm cap
White

CP2078/01

Connects different parts of your breast pump

Find matching products on the specifications tab

This slave cap links the breast pump body to the motor unit. It is suitable for the electric pump. Package comes with one cap.

PHILIPS

Specifications

Replaceable parts

- Fit product types: SCF301/01 SCF301/02 SCF301/

03 SCF301/04 SCF301/05 SCF303/01 SCF303/02



Issue date 2023-03-08

Version: 1.1.1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com