



Philips Zoom Chairside
In office No Touch Up -2pk

With Relief ACP Oral Care Gel

DIS570/01

Philips Zoom In-Office Procedure Kits

For use with WhiteSpeed Lamp and Procedure gels

With the input of hundreds of Philips Zoom users, Philips has made administering the Philips Zoom whitening system easier, faster and safer than ever before.

Your most confident smile now

- Light-activated technology for powerful whitening

Impressive results fast

- White bright smile to improve your patient's confidence

Includes ACP Oral Care gel

- Relief[®] ACP Oral Care Gel helps manage patient's sensitivity

PHILIPS ZOOM!

Specifications

Consumer Trade Item

- Height: 6
- Width: 17.5
- Length: 26.2
- Net Weight: 0.239
- Gross Weight: 0.36
- GTIN: 08710103839361
- Country of origin: US
- Harmonized system code: 33069000

Outer Carton

- Height: 34.9
- Width: 41.3
- Length: 46.4
- Net Weight: 2.868
- Gross Weight: 5.068
- GTIN: 18710103839368

Performance

- Whitening results: Maximum whiteness in 45 min.

Highlights

Includes ACP Oral Care gel

Recommendation is to place a light coating of Relief ACP (amorphous calcium phosphate) Oral Care Gel on the lingual surfaces of teeth on both arches before whitening treatment and to apply ACP immediately following an in-chair procedure directly onto the teeth or in custom trays for 10 to 30 minutes

Impressive results fast

Your patients can get up to eight shades whiter teeth in 45 minutes, excluding preparation time

Light-activated technology

Unlike many other systems, Philips Zoom WhiteSpeed uses a proprietary LED light to accelerate whitening for maximum results in minimal time.



Issue date 2023-03-16

Version: 7.7.1

EAN: 87 10103 83936 1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com