

PHILIPS

AVENT



Electric breast pump



SCF395

Inspired by baby, made for mum

Enter a new era of expression with the Philips Avent Electric Breast Pump. Perfect balance of suction & nipple stimulation, inspired by baby's drinking rhythm and soft cushion gently adapts to your nipple size & shape. See more below.

Other benefits

- Choose the right teat flow for your baby

Personalised

- Minimal parts and intuitive setup
- Quiet motor for a discreet experience anywhere (corded use)
- Start/Pause function
- Personalised experience, 8 + 16 setting levels
- Memory function

Gentle and comfortable

- Soft and adaptive silicone cushion
- Express without leaning forward


Quick expressing, more milk in less time*

- Natural motion technology for a quick milk flow*

Finding the right teat is important


Highlights

Quick milk flow*



Express more milk in less time* with a cushion that stimulates the breast to express milk just like baby. It seamlessly adjusts from stimulation mode to expression mode and applies just the right amount of teat stimulation and suction for maximum milk flow. Based on milk flow initiation time (time to Milk Ejection Reflex — MER) results.*

Adapts to your nipple



Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your nipple. It fits 99.98% of nipple sizes* (up to 30 mm).

Finding the right teat is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a teat with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

Sit comfortably



Feel relaxed as you express thanks to a design that lets you sit upright rather than having to lean forwards. No need to lean forwards — clinically proven comfortable expressing position*. Based on product clinical trial results with 20 participants (2019); 90% of participants find the expressing position comfortable (single electric); 95% of participants find the expressing position comfortable (double electric).

Personalise to your needs



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 8 stimulation and 16 expression levels for a personalised experience.

Express anywhere discreetly



Express whenever and wherever you like thanks to our quiet motor with a slimline design. Corded use only.

Easy to clean and set up



Our closed expression system means milk stays out of the tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Highlights

Pause whenever you like



If you need to adjust or take a break, the start/pause button is always at your fingertips.

Save your favourite settings



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Find the right teat flow



Every baby feeds differently and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalise your bottle. All of the Natural Response Teats are made of soft silicone.

Specifications

Power

Mains voltage: 100 - 240 V

What is included

Expression kit: 1 pcs

Motor unit (Corded use only): 1 pcs

Sealing disc: 1 pcs

Micro-USB adapter (corded use): 1 pcs

Milk bags: 5 pcs

125 ml (4 oz) Bottle: 1 pcs

Slow Teat: 1 pcs

