

# PHILIPS

## sonicare

Two sonic electric  
toothbrushes

Essence

1 mode

2 brush heads

HX5352/46



## Electric toothbrush for better plaque removal

The unique dynamic cleaning action of Philips Sonicare HX5352/46 gently and effectively reaches deep between teeth and along the gumline.

### **Proven to improve oral health**

- Clinically proven safe and gentle

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time

### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.




## SmarTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: NiMH

## Design and finishing

Color: White

## Service

Warranty: 2-year limited warranty

## Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy screw-on brush heads

Handle: Ergonomic design, Rubber grip for easy handling

Brushing time: Up to 2 weeks

## Items included

Charger: 1

Handles: 2 Essence

Travel case: 1

Brush heads: 2 e-Series standard

## Cleaning performance

Performance: Up to 2X more plaque removal\*

Health benefits: Helps improve gum health

Speed: Up to 31000 brush movement/min

Timer: SmarTimer

Whitening benefits: Helps naturally remove stains

## Modes

Clean: For exceptional everyday clean

