

Philips  
Waist massager

### Combined Massage

Stereo support  
Wireless use

PPM4721



## Multifunctional massage, release pressure

Combining 3 functions to energize the waist

Made of light and thin fabric, this product can relieve pain through vibration and hot compress. The 3D support bar can stabilize the waist to prevent excessive bending. Equipped with a cold gel bag for calming and relieving muscle sprains.

#### High-frequency vibration, relieve waist discomfort

- 4-zone high-frequency vibration massage to relax and relieve
- 4 vibration modes to enrich massage experience

#### 3D strong support, reduce lumbar pressure

- Metal steel bar support bar to disperse support pressure

#### Smart massage, safe to use

- 40-55°C hot compress temperature, 4 adjustable temperatures
- 30-minute intelligent timed massage, avoid excessive massage
- Built-in battery, avoid power cord constraints.
- Soft touch, light fabric

# PHILIPS

Waist massager  
Combined Massage Stereo support, Wireless use

# Specifications

## General specifications

- Color: Dark blue
- Net weight: approx. 1.1kg
- Product Dimension (LxWxH): 1090x224x50 mm
- Main materials: Mesh cloth

## Main parameters

- Working voltage: 3.7 V
- Temperature of heating: 40-55 °C
- Auto shut-off time: 30 minutes
- Working power: 9.5 W

- Battery: 3200 mAh
- Charging current: 2 A
- Charging voltage: 5 V
- Charging port: USB Type-C

## Country of origin

- Made in: China

## Service

- Warranty period: 2 years

PPM4721DB/97

# Highlights

## High frequency vibration

Built-in 4 high-frequency vibration massage areas, directly hit the sore muscles. You can relax your waist muscles immersively and relieve muscle fatigue and soreness.

## Multiple vibration modes

There are 2 warm-up massage modes and 2 recovery massage modes after exercise. With 4 zones of cool rhythmic lighting, the massage has a more technological sense.

## 3D strong support function

There are 2 arch-shaped supporting steel bars to stabilize the lumbar spine, disperse the pressure on the lumbar spine from the upper body, and provide protection and fixation to the waist.

## Hot compress function

NTC intelligent temperature control can strengthen the local blood circulation of the human body, thereby effectively relieving muscle tension and pain, allowing you to relax and relieve the fatigue of the day.

## Timing shut-off

A timed 30-minute working time for automatically shut off, effectively avoiding excessive massage and causing discomfort of redness and swelling of the skin.

## Use wireless

The control box has a built-in battery with a capacity of 3200mAh, which can be used for 120 minutes when fully charged.

## Skin-friendly and soft fabric

It is made of fine fabric sewing, and the whole machine feels soft to the touch. The fabric is added with antibacterial and antibacterial technology, which makes it healthier to wear and more secure to use.



Issue date 2023-05-15

Version: 1.1.1

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)