

Philips Sonicare  
Rechargeable powered  
dental flosser - silver

**2 silver nozzles**  
237 ml mouthwash



**HX8472/11**

## Healthier gums in 2 weeks, guaranteed\*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro dental flosser is the easiest way to effectively clean between teeth. AirFloss Pro can be used with mouthwash or water and is clinically proven to be as effective as floss for gum health.\*\*

### Improves oral health

- Up to 99.9% plaque removal\*\*\*
- Helps prevent cavities between teeth
- Clinically proven to be as effective as floss for gum health\*\*

### Superior technology

- Air and micro-droplet technology
- Triple-burst customisable settings
- High-performance nozzle

### The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- Total daily confidence with a fresh and clean mouth
- An easy way to start a healthy habit

**PHILIPS**  
**sonicare**


# Highlights

## Better gum health



Philips Sonicare AirFloss Pro/Ultra nozzle is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

## Micro-droplet technology



Our clinically proven results are possible thanks to our unique technology, which combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## High-performance nozzle



The new AirFloss Pro/Ultra high-performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple) and hold down the activation


button for continuous automatic bursts or press and release for manual burst mode.

## Fresh and clean with mouthwash




Fill the reservoir on the handle with either mouthwash or water, and then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Triple-burst customisation



AirFloss Pro/Ultra's powerful spray can be customised to your liking, with single, double or triple bursts for each press of the activation button.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Consumer Trade Item

- Height: 23.5 cm
- Width: 8.4 cm
- Length: 15.7 cm
- Net weight: 0.59 kg
- Gross weight: 0.714 kg
- GTIN: 8710103817079
- Country of origin: RO
- Harmonized system code: 85098000

## Outer Carton

- Height: 25.5 cm
- Width: 16.9 cm
- Length: 27 cm
- Net weight: 1.77 kg
- Gross weight: 2.363 kg
- GTIN: 18710103817076

## Design and finishing

- Colour: White with grey accents

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Battery Life: 1–3 weeks\*\*\*\*
- Nozzle attachment: Easily snaps on and off

## Items included

- Charger: 1
- Mouthwash: 237 ml BreathRx
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 2

## Cleaning performance

- Health benefits: Cleans entire mouth in 60 secs, Improve gum health in 2 weeks\*
- Performance: Removes up to 99.9% plaque\*\*
- Clinically proven results: As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction
- Customisable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
- For best results: Change nozzle every 6 months

## Ingredients

- Mouth Rinse Ingredients: Aqua, Sorbitol, Propylene Glycol, PEG-40, Hydrogenated Castor Oil, Xylitol, Poloxamer 407, Aroma, Zinc Gluconate, Cocamidopropyl Betaine, Cetylpyridium Chloride, Sodium Saccharin, Citric Acid, Limonene, BHT, CI42090



Issue date 2024-09-02

Version: 2.0.1

EAN: 87 10103 81707 9

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* or your money back

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouthwash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under the Support tab for further details.

\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary

\*\*\*\* depending on burst setting used

\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.