

Philips  
Massage Headband

6-area kneading massage

Colorful screen controller  
Adjustable headband  
Breathable fabric

PPM1321



## Enjoy a relaxing head massage

With kneading and hot compress

This band-type head massager provides a new way of head massage. It can comfortably fit different head sizes. With its portable size, you can enjoy relaxing kneading head massage and hot compress anywhere.

### Comfortable head massage experience

- 6-area air pressure kneading massage
- 3 levels of constant warm compress

### Designed for everyday massage

- Portable design
- Adjustable design to fit different head circumference
- Suitable for use before sleep

### Easy to control

- Simple user interface

# PHILIPS

## Massage Headband

6-area kneading massage Colorful screen controller, Adjustable headband, Breathable fabric

PPM1321LB/97

# Specifications

### General specifications

- Color: Light Blue
- Net weight: 510 g
- Massage type: Kneading, Hot compress
- Packaging Dimension(L)\*(W)\*(H): 346x226x48 mm

### Country of origin

- Country of origin: China

### Service

- Warranty period: 2 years

# Highlights

### 6-area air pressure kneading massage

The massage headband fits the head circumference and has 6 massage areas covering the forehead, temples, and occipital region. The air bag simulates kneading pressure of human hand massage, which can effectively relax the head muscles and enhance blood circulation.

### 3 levels of constant warm compress

The hot compress module with built-in temperature sensor can provide hot compress at a well-controlled temperature range, which can help to improve the blood circulation at forehead area.

### Portable design

Unlike other head massager with helmet shape and size, this headband device is super light weight and portable. The combination of fabric main unit with external controller makes it easy to be folded to a handy size and take away.

### Adjustable to fit different head circumference

With the Velcro at back of the headband. It can fit properly with most head circumference of adults. Which is essential for best head massage experience. Users can also adjust it easily when they feel it too tight or too loose.

### Suitable for use before sleep

When you enjoy a comfortable massage with this device and gradually fall asleep on bed, the lightweight headband will not bring any discomfort.

### Simple user interface

There are only 2 buttons, one for massage modes, another for hot compress levels. The color screen around the buttons directly shows the status. No complex process, just click and enjoy it.



Issue date 2023-05-22

Version: 1.1.1

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)