

# PHILIPS

# AVENT



Baby Bottle with  
Airfree vent

Natural Response

2 Bottles

4oz/125ml

Slow Flow Nipple

0-3m



SCY670/02

## Flows like a breast, aids to reduce feeding issues

### A nipple that works like a breast\*

Provides calm, comfortable feeds. The Natural Response Teat supports baby's unique drinking rhythm, while the AirFree vent is designed for extra protection against colic, gas and reflux. Finding the right nipple is important. See more below.

#### Easy to combine breast and bottle feeding

- Nipple releases milk when baby actively drinks
- Natural latch on with breast-shaped nipple

#### Other benefits

- Simple to use, easy to clean and quick to assemble
- Easy to hold even for little hands
- Choose the right nipple flow for your baby
- Same products, new navigation
- Explore our bottle range
- No-drip nipple design prevents spills and lost milk
- Compatible across the Philips Avent range
- Natural Response Nipples and Bottles are BPA free\*
- Be patient as baby adjusts

#### AirFree vent designed for less air ingestion

#### Finding the right nipple is important

# Natural Response

Baby Bottle with Airfree vent

## Highlights

### Natural latch on



The wide, soft and flexible nipple is designed to mimic the shape and feel of a breast, helping baby to latch on and feed comfortably.

### Easy to hold



The ergonomic bottle is easy to grip at any angle to give maximum comfort during feeding. Easy for your hands and tiny hands to hold.

### Finding the right nipple is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a nipple with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

### Simple to use and clean




Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly, as AirFree vent is a single piece.

### Natural Response Nipple



The Natural Response Nipple works with your baby's natural feeding rhythm, making it easy to combine breast feeding and bottle feeding. The nipple has a unique opening which only releases milk when the baby actively drinks. So when they pause to swallow and breathe, the milk pauses too.

### Airfree vent



The AirFree vent is designed to offer extra protection against feeding issues by preventing air from getting into your little ones tummy while upright feeding.

SCY670/02

### Find the right nipple flow



Every baby feeds differently, and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalize your bottle. All of the Natural Response Nipples are made of soft silicone.

### Same products, new navigation



We have moved to a pace-based flow navigation system. Start with the nipple that comes with the bottle. Try a lower flow if milk is leaking from baby's mouth or baby is gulping. Try a higher flow if baby is playing with the nipple instead of drinking or seems frustrated. As we make this update, you may receive either style of pack.

### No-drip nipple



The nipple opening is designed to release milk only when the baby is feeding. So you can confidently avoid milk loss whether at home or on the go.

# Natural Response

Baby Bottle with Airfree vent

SCY670/02

## Highlights

### Compatible across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it.

**BPA free\***



The Philips Avent Natural bottles and nipples are made of BPA free\* material.

### Be patient as baby adjusts



Our new Natural Response Baby Bottles are different from free-flow baby bottles. Just like breastfeeding, it could take a few tries to get it right. That's perfectly natural.

### Explore our bottle range



Our bottles are available in different sizes to better suit your baby's needs. Choose from a range of options to find the perfect fit for your little one, ensuring comfort and convenience at every stage.

# Natural Response

Baby Bottle with Airfree vent

SCY670/02

## Specifications

### Design

**Bottle design:** Transparent

### Material

**Bottle:** Polypropylene, BPA free\*

**Nipple:** Silicone, BPA free\*

### What is included

**AirFree vent:** 2 pcs

**4oz/125ml Baby Bottle:** 2 pcs

**Slow Nipple:** 2 pcs

### Functions

**Nipple Features:** Natural latch on, No-drip design, Soft and flexible

**Bottle design:** Wide neck, Ergonomic shape

**Bottle ease of use:** Easy to assemble, Easy to clean, Dishwasher safe, Easy to hold

### Development stages

**Stage:** 0-3m

