



PHILIPS

Travel kit for DreamStation
CPAP device

DreamStation

Case color: light gray
Including 65W Power Supply
with USB port

 Product Picture

HH1423

Ideal for travelling

Compact travel bag for DreamStation

Gives patients a convenient, compact, and discreet way to travel with DreamStation

Ideal combination of protection and portability

- Benefit from your CPAP therapy on the move

Ideal for travel

- Dedicated compartments

Including 65W Power Supply with USB Port

- with USB Port

Highlights

Lightweight power supply



The case includes a 65W power supply which is smaller than the power supply used for the DreamStation machine and humidifier. The power supply has a built-in USB port for charging cell phones and other electronics while using the DreamStation machine.

Compact and easy to use

A molded insert liner helps the patient thoroughly pack the PAP system components: mask and tubing, 65W power supply, DreamStation WITHOUT humidifier.

Easy to care and sturdy

Lightweight, hard-fabric case. Removable shoulder strap + trolley back strap to easily slide onto luggage handles. Secure zippered closures help facilitate airport screening; and the travel case easily fits in airplane overhead bins and under airline seats.

Specifications

Travel

Dedicated to DreamStation: Travel case with 65W power supply with USB port holds mask, headgear, tubing and device

Service

Warranty: 90-days

