

Avent Breast Pumps
Silicone diaphragm

Diaphragm

CP2084/01



To create a vacuum in the breast pump.

Find matching products on the specifications tab

This silicone diaphragm creates a vacuum in your manual breast pump.

AVENT

Specifications

Replaceable parts

- Fit product types: SCF323, SCF363, SCF903/01



Issue date 2023-07-15

Version: 2.2.1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com