

# PHILIPS

Under-the-nose  
mask

Nuance Gel

Tube-in-front design

Gel pillows cushion

Choose between fabric or gel  
padded frame

Available in multiple cushion  
sizes



HH1100/00



## Nasal gel pillow mask

The Nuance gel mask platform allows you to choose between a fabric frame or gel padded frame for optimal comfort.

### Easy to use

- Light-weight, flexible tubing for less pull on the mask.

### Greater stability

- Stays in place throughout the night.

# Highlights

## Lightweight and simple



Nuance Gel has a light-weight, flexible tubing that allows for easy assembly and disassembly.

## Innovative non-slip headgear



Ideal for patients that want to sleep freely.

# Specifications

## Specifications

Operating pressure: 4 -20 cmH2O

Sound pressure level: 19dBA

N° of parts: Four

Cushion material: Silicone rubber

Headgear material: Nylon UBL/Urethane foam/Lycra

Frame material: Nylon, Lycra, Urethane foam and Trico mesh

Cleaning instructions: **mask:** Hand wash with mild detergent, Rinse thoroughly, air dry

**Cleaning: headgear, tubing:** Hand wash with mild detergent

**Swivel/Tubing/Elbow material:** Polypropylene

**Pillows material:** Urethane gel

## Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

## Design

**Easy to use:** Easy assembly and disassembly

**Comfort:** Fabric frame

## Service

**Warranty:** 90-days

