

# PHILIPS

Mask - Fit Pack (S,  
M, L)

Nuance Pro

Gel

Pack

3 pillows included (S, M, L)



HH1101



## A new choice for better sleep

The gel nasal pillow technology of Nuance delivers exceptional comfort and seal, helping you to support better sleep.

### Easy to use

- Light-weight, flexible tubing for less pull on the mask.

### Gel pillows

- The pillows fit difference nostril sizes.

### Greater stability

- Stays in place throughout the night.

# Highlights

## Lightweight and simple



Easier to assemble, disassemble and use than the leading pillows mask.

## Provides comfort and seal



Nuance gel pillows cushion is designed differently to provide a more comfortable experience. The gel base of the pillows should fit comfortably under your nose.

## Innovative non-slip headgear



Ideal for patients that want to sleep freely.

# Specifications

## Specifications

**Operating pressure:** 4 -20 cmH<sub>2</sub>O

**Sound pressure level:** 19 dBA

**No. of parts:** Four

**Cushion material:** Silicone rubber

**Pillows/frame pads material:** Urethane gel

**Headgear material:** Nylon UBL/Urethane foam/Lycra

**Frame material:** Thermoplastic polyester elastomer

**Cleaning instructions: mask:** Hand wash with mild detergent, Rinse thoroughly, air dry

**Cleaning: headgear, tubing:** Hand wash with mild detergent

**Swivel/Tubing/Elbow material:** Polypropylene

## Intent of use

To be used on patients <gt/>30 kg: with CPAP prescription

## Design

**Easy to use:** Easy assembly and disassembly

## Service

**Warranty:** 90-days

