

# PHILIPS

Sonic electric  
toothbrush with app

DiamondClean  
9000

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities

HX9911/75



## Whiter, healthier teeth for life

Whiter teeth in just 1 day\*\*\*\*

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

### Enjoy a stylish and user-friendly design

- Simple to charge and style to boot

### DiamondClean users experience improved oral health

- Deepest clean while protecting your sensitive gums
- Up to 7x healthier gums\* in just 2 weeks

### Let the DiamondClean 9000 adapt to your needs

- Personalized brushing experience
- Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- Always know when to replace your brush heads

### Stay on track with a personalized progress report

- Start and keep up healthy habits

# Highlights

## Up to 7x healthier gums\*




Focus on improving your gum health with the G3 Premium Gum Care brush head. Its smaller size and targeted gum line bristles make it the perfect way to keep this area gum disease free. People who use this brush head have up to 100% less gum inflammation\* and up to 7x healthier gums in just two weeks.\*

## G3 Premium Gum Care brush head



The G3 Premium Gum Care brush head is designed to help you improve your gum health. Its smaller size and targeted gumline bristles gently and effectively clean the gum line, where gum disease commonly starts out.

## Progress Report



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the

Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

## 4 modes, 3 intensities




The DiamondClean 9000 comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

## BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

## Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

## BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

## Long-life battery



Your DiamondClean 9000 will last for up to two weeks on a single charge with regular use.

# Specifications

## Connectivity

**Bluetooth® wireless technology:** Connected brushing app

## Power

**Voltage:** 110–220 V

## Technical specifications

**Operating time (full to empty):** 14 days\*\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

## Design and finishing

**Color:** Black

## Service

**Warranty:** 2-year limited warranty

## Compatibility

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

## Ease of use

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Replacement reminder:** To always ensure best results, reminder icon lights up

## Items included

**Brush heads:** 1 G3 Premium Gum Care

**Travel case:** 1

**Handle:** 1 DiamondClean 9000

**Glass charger:** 1

**Charger base:** 1

## Cleaning performance

**Health benefits:** Up to 7x healthier gums\*

**Whitening benefits:** Whitens teeth in 1 day\*\*\*\*

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Vibrates handle to alert user

## Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**White+:** To remove surface stains

## Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

**BrushSync Technology:** Connects smart handle and, smart brush head



\* more than a manual toothbrush

\*\* compared to DiamondClean

\*\*\*\* in White+ Mode with leading whitening toothpaste vs a manual toothbrush

\*\*\* based on two periods of two-minute brushings per day