

PHILIPS

AVENT

Essential

Single Electric
Breast Pump

SCF323/11



Express effortlessly

This Philips Avent single electric breast pump keeps milk flow at an optimum and gently adapts to your size and shape for a more comfortable experience. It also comes equipped with the latest USB-C format for fast charging.

Personalized

- Personalised experience, 4 + 4 setting levels
- Minimal parts and intuitive setup
- Memory function
- Small and Compact

Gentle and comfortable

- Soft and adaptive silicone cushion


Highlights

Soft adaptive silicone cushion



One size fits all. The silicone cushion gently flexes and adapts to fit your nipple.

Personalised experience



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 4 stimulation and 4 expression levels for a personalised experience.

Easy to clean and set up



Our closed expression system means milk stays out of the tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Memory setting



Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Small and Compact



Small and Compact motor unit

Specifications

Power

USB-C connection

This appliance does not contain a power adapter.: Yes. Only USB-C cable provided.

Main voltage: 5V and ≥ 1 A

Supply unit is not included. Not chargeable, corded use only.

Material

Bottle: BPA free*

Teat: Silicone

0% BPA*

What is included

Motor unit

Silicone tube

USB-C cable

Natural bottle (125 ml)

Expression kit

Screw-on lid

Ease of use

Breast pump use: Easy cleaning & assembly, memory function, rechargeable battery

Functions

Settings: 4 expression levels, 4 massage levels

