



HeartStart Intrepid Monitor/Defibrillator

P/N: 867172

This easy-to-use, lightweight and rugged monitor/defibrillator gives you the power to respond quickly and act confidently when your patient's life is on the line. Other relevant proof points are included in the five focal areas below.



Our EcoPassport

As a company committed to doing business sustainably, we are keen to help our customers make responsible choices. We offer solutions that improve people's health and well-being while reducing impact on the environment.

Our EcoPassports summarize the environmental benefits our products offer in one or more of our focal areas. For example, increased energy efficiency, more sustainable packaging, or a circular-ready product design; optimized for repair, refurbishing and recycling.

In this way, we want to help ensure that each purchase decision is the right one for our customer's needs *and* the planet.

Further reading at:
www.philips.com/sustainability



Energy

- Power in on / scan mode: 52.5 W
- Power for battery charging: 46.2 W
- Power in standby-mode: 0.076 W
- Annual energy consumption: 157.63 kWh***



Packaging

- **Total Weight: 2.4 – 2.7 kg**
- Cardboard: 2.25-2.55kg
- LDPE bag: 0.15 kg
- Recycled cardboard: 50%
- No use of polyvinyl chloride (PVC) or expanded polystyrene (EPS)



Weight

- Product: 6.9 kg*



Substances

- RoHS 2 compliant**



Circularity

- Service and spare parts available

* Base unit includes pads adaptor cable, full roll paper, battery pack, ac power supply

**Directive 2011/65/EU plus amendment 2015/863

*** Average use scenario